










Prima Plus™ Implant Drill Usage Chart

The drills for Prima Plus™ Implants utilize an improved design for increased efficiency, longevity and overall ease-of-use. To maintain optimal performance Keystone Dental recommends replacing drills after 20 uses (depending upon bone density). The chart below provides an easy method for tracking drill usage. Simply check the box next to the corresponding drill each time it is used. Then, order a replacement drill as you near 20 uses.

						REORDER	CATALOG #
FINAL TWIST DRILLS		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15750K
	1.8mm Initial Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15751K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15752K
	2.2/2.6mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15753K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15754K
	2.6/2.9mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15755K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15756K
	2.9/3.3mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15757K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15758K
	3.4/3.9mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15759K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15760K
	3.8/4.1mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15761K
		8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15762K
	4.4/4.8mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15763K
	8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15764K	
4.6/5.3mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15765K	
	8-15mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15766K	
5.3/5.7mm Twist Step Drill	8-11.5mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	15767K	